

## SPACA3 protein (STJ505365)

STJ505365

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	SPACA3-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">124912</a>
<b>Gene Symbol</b>	<a href="#">SPACA3</a>
<b>Uniprot ID</b>	<a href="#">SACA3_HUMAN</a>
<b>Immunogen</b>	Synthetic peptides taken within amino acid region 80-130 on human sperm acrosome membrane-associated protein 3 isoform 2.
<b>Immunogen Region</b>	80-130
<b>Specificity</b>	
<b>Immunogen Sequence</b>	