

GTPase HRAS protein (STJ505370)

STJ505370

GENERAL INFORMATION

Product Type	Proteins
Short Description	GTPase HRAS-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3265
Gene Symbol	HRAS
Uniprot ID	RASH_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 139-189 on human GTPase HRAS protein.
Immunogen Region	139-189
Specificity	
Immunogen Sequence	