

PAX2 isoform a + e protein (STJ505371)

STJ505371

GENERAL INFORMATION

Product Type	Proteins
Short Description	PAX2 isoform a + e-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5076
Gene Symbol	PAX2
Uniprot ID	PAX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-220 on human PAX2 isoform d protein.
Immunogen Region	200-220
Specificity	
Immunogen Sequence	