

**cKit protein (STJ505376)**  
STJ505376**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | cKit-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | 3815   |
| <b>Gene Symbol</b>        | KIT  |
| <b>Uniprot ID</b>         | KIT_HUMAN  |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 915-965 and 950-976 on human Tyrosine-Kinase c-Kit protein. |
| <b>Immunogen Region</b>   | 915-965 and 950-976  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |