

hCG beta Peptide protein (STJ505378)

STJ505378

GENERAL INFORMATION

Product Type	Proteins
Short Description	hCG beta Peptide-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114335
Gene Symbol	CGB1
Uniprot ID	CGB1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 142-165 on human Chorionic gonadotropin beta protein.
Immunogen Region	142-165
Specificity	
Immunogen Sequence	