

hCG beta 2 Peptide protein (STJ505380)

STJ505380

GENERAL INFORMATION

Product Type	Proteins
Short Description	hCG beta 2 Peptide-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 82-132 on human choriogonadotropin subunit beta variant 2 protein.
Immunogen Region	82-132
Specificity	
Immunogen Sequence	