

ErbB 2 protein (STJ505383)

STJ505383

GENERAL INFORMATION

Product Type	Proteins
Short Description	ErbB 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2064
Gene Symbol	ERBB2
Uniprot ID	ERBB2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1210-1255 receptor tyrosine-protein kinase erbB-2 isoform a.
Immunogen Region	1210-1255
Specificity	
Immunogen Sequence	