

Zika Virus NS5 protein (STJ505384)

STJ505384

GENERAL INFORMATION

Product Type	Proteins
Short Description	Zika Virus NS5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic pepitde taken within amino acid region 3230-3280 on Zika virus NS5 protein.
Immunogen Region	3230-3280
Specificity	
Immunogen Sequence	