

KK-LC-1 protein (STJ505392)

STJ505392

GENERAL INFORMATION

Product Type	Proteins
Short Description	KK-LC-1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	203413
Gene Symbol	CT83
Uniprot ID	KKLC1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 20-70 and 61-111 on human Kita-kyushu lung cancer antigen 1.
Immunogen Region	20-70 and 61-111
Specificity	
Immunogen Sequence	