

YME1L1 protein (STJ505393)

STJ505393

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | YME1L1-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 10730 |
| Gene Symbol | YME1L1 |
| Uniprot ID | YMEL1_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 720-773 on human ATP-dependent metalloprotease YME1L1 protein. |
| Immunogen Region | 720-773 |
| Specificity | |
| Immunogen Sequence | |