

MYLK protein (STJ505395)

STJ505395

GENERAL INFORMATION

Product Type	Proteins
Short Description	MYLK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4638
Gene Symbol	MYLK
Uniprot ID	MYLK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1870-1917 on human Myosin light chain kinase, smooth muscle isoform 1 protein.
Immunogen Region	1870-1917
Specificity	
Immunogen Sequence	