

MOR1 protein (STJ505397)
STJ505397**GENERAL INFORMATION**

Product Type	Proteins
Short Description	MOR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4988
Gene Symbol	OPRM1
Uniprot ID	OPRM_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 405-455 on human Mu-type opioid receptor isoform MOR-1i protein.
Immunogen Region	405-455
Specificity	
Immunogen Sequence	