

DWORF protein (STJ505399)

STJ505399

GENERAL INFORMATION

Product Type	Proteins
Short Description	DWORF-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Strit1
Uniprot ID	DWORF_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-30 on mouse Sarcoplasmic/endoplasmic reticulum calcium ATPase regulator DWORF protein.
Immunogen Region	1-30
Specificity	
Immunogen Sequence	