

B-raf protein (STJ505403)

STJ505403

GENERAL INFORMATION

Product Type	Proteins
Short Description	B-raf-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	673
Gene Symbol	BRAF
Uniprot ID	BRAF_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-750 on human serine/threonine-protein kinase B-raf protein.
Immunogen Region	700-750
Specificity	
Immunogen Sequence	