

GPCR GPR8 protein (STJ505407)

STJ505407

GENERAL INFORMATION

| | |
|--------------------------|--------------------------------------------------------|
| Product Type | Proteins |
| Short Description | GPCR GPR8-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|-----------------------------------------------------------------|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: IHC: 1:50-1:100 |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------------------------------------------------------------------------------------|
| Gene ID | 2832 |
| Gene Symbol | NPBWR2 |
| Uniprot ID | NPBW2_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 280-333 on human Neuropeptides B/W receptor type 2 protein. |
| Immunogen Region | 280-333 |
| Specificity | |
| Immunogen Sequence | |