

STAT4 protein (Phospho) (STJ505411)

STJ505411

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-STAT4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6775
Gene Symbol	STAT4
Uniprot ID	STAT4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 690-740 on human Signal transducer and activator of transcription 4 protein.
Immunogen Region	Phosphorylated serine 721 residue.
Immunogen Specificity	690-740
Immunogen Sequence	