

STAT6 protein (STJ505413)
STJ505413**GENERAL INFORMATION**

Product Type	Proteins
Short Description	STAT6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6778
Gene Symbol	STAT6
Uniprot ID	STAT6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 655-705 on human Signal transducer and activator of transcription 6 protein.
Immunogen Region	655-705
Specificity	
Immunogen	
Sequence	