

## JAK2 protein (Phospho) (STJ505417)

STJ505417

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Phospho-JAK2-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3717</a>
<b>Gene Symbol</b>	<a href="#">JAK2</a>
<b>Uniprot ID</b>	<a href="#">JAK2_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 980-1030 on human JAK2 protein. Phosphorylated 1008 & 1009 tyrosine residues.
<b>Immunogen Region</b>	980-1030
<b>Specificity</b>	
<b>Immunogen Sequence</b>	