

JAK3 protein (Phospho) (STJ505419)

STJ505419

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-JAK3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3718
Gene Symbol	JAK3
Uniprot ID	JAK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 750-800 on human JAK3 protein. Phosphorylated tyrosine 785 residue.
Immunogen Region	750-800
Specificity	
Immunogen Sequence	