

MAPK11 protein (STJ505433)

STJ505433

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAPK11-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:50-1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5600
Gene Symbol	MAPK11
Uniprot ID	MK11_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Mitogen-activated protein kinase 11.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	