

MAPK2 protein (STJ505440)

STJ505440

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAPK2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9261
Gene Symbol	MAPKAPK2
Uniprot ID	MAPK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Mitogen-activated protein kinase 2.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	