

MAPK6 protein (STJ505446)

STJ505446

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAPK6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5597
Gene Symbol	MAPK6
Uniprot ID	MK06_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 480-530 on human Mitogen-activated protein kinase 6.
Immunogen Region	480-530
Specificity	
Immunogen Sequence	