

MAPK7 protein (STJ505452)

STJ505452

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAPK7-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5598
Gene Symbol	MAPK7
Uniprot ID	MK07_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-750 on human Mitogen-activated protein kinase 7.
Immunogen Region	700-750
Specificity	
Immunogen Sequence	