

Zika Virus Envelope Protein E protein (STJ505454)

STJ505454

GENERAL INFORMATION

Product Type	Proteins
Short Description	Zika Virus Envelope Protein E-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 280-320 on Zika Virus Envelope Protein E.
Immunogen Region	280-320
Specificity	
Immunogen Sequence	