

Zika Virus NS1 protein (STJ505463)

STJ505463

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | Zika Virus NS1-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:250-1:500 |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |