

MAVS protein (STJ505464)

STJ505464

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAVS-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57506
Gene Symbol	MAVS
Uniprot ID	MAVS_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-540 on on human mitochondrial antiviral-signaling protein isoform 1.
Immunogen Region	500-540
Specificity	
Immunogen Sequence	