

PD-L1 protein (STJ505468)
STJ505468**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PD-L1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
Uniprot ID	PD1L1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 80-130 on human Programmed cell death 1 ligand 1 protein.
Immunogen Region	80-130
Specificity	
Immunogen Sequence	