

PD-L2 protein (STJ505470)
STJ505470**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PD-L2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80380
Gene Symbol	PDCD1LG2
Uniprot ID	PD1L2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 115-165 on human Programmed cell death 1 ligand 2 protein.
Immunogen Region	115-165
Specificity	
Immunogen Sequence	