

ABHD2 protein (STJ505472)
STJ505472**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ABHD2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54608
Gene Symbol	Abhd2
Uniprot ID	ABHD2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse a/b hydrolase II protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	