

ZFP36L1 protein (STJ505474)

STJ505474

GENERAL INFORMATION

Product Type	Proteins
Short Description	ZFP36L1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	677
Gene Symbol	ZFP36L1
Uniprot ID	TISB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Zinc finger protein 36, C3H1 type-like 1.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	