

## ZFP36L2 protein (STJ505481)

STJ505481

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Proteins   |
| <b>Short Description</b> | ZFP36L2-protein if for use in research applications. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       |  |
| <b>Reactivity</b>        |  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | WB: 1:500   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 678        |
| <b>Gene Symbol</b>        | ZFP36L2    |
| <b>Uniprot ID</b>         | TISD_HUMAN |
| <b>Immunogen</b>          |            |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |