

IL18RAP protein (STJ505482)

STJ505482

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL18RAP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8807
Gene Symbol	IL18RAP
Uniprot ID	I18RA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Interleukin-18 receptor accessory protein.
Immunogen Region	500-550
Specificity	
Immunogen Sequence	