

**CX3CL1 protein (STJ505484)**  
STJ505484**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	CX3CL1-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">20312</a>
<b>Gene Symbol</b>	<a href="#">Cx3cl1</a>
<b>Uniprot ID</b>	<a href="#">X3CL1_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 350-395 on mouse Fractalkine protein.
<b>Immunogen Region</b>	350-395
<b>Specificity</b>	
<b>Immunogen Sequence</b>	