

BMPR1A protein (STJ505490)

STJ505490

GENERAL INFORMATION

Product Type	Proteins
Short Description	BMPR1A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:50-1:150
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	657
Gene Symbol	BMPR1A
Uniprot ID	BMR1A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 460-510 on human Bone morphogenetic protein receptor type-1A.
Immunogen Region	460-510
Specificity	
Immunogen Sequence	