

IQSEC3 protein (STJ505492)

STJ505492

GENERAL INFORMATION

Product Type	Proteins
Short Description	IQSEC3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	440073
Gene Symbol	IQSEC3
Uniprot ID	IQEC3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1100-1150 on human IQ motif and SEC7 domain-containing protein 3.
Immunogen Region	1100-1150
Specificity	
Immunogen Sequence	