

SQSTM1/p62 protein (STJ505496)

STJ505496

GENERAL INFORMATION

Product Type	Proteins
Short Description	SQSTM1/p62-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8878
Gene Symbol	SQSTM1
Uniprot ID	SQSTM_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-440 on human Sequestosome-1 protein.
Immunogen Region	400-440
Specificity	
Immunogen Sequence	