

PRAS40 protein (STJ505498)
STJ505498**GENERAL INFORMATION**

Product Type	Proteins
Short Description	PRAS40-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84335
Gene Symbol	AKT1S1
Uniprot ID	AKTS1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 160-210 on human Proline-rich AKT1 substrate 1 isoform a protein.
Immunogen Region	160-210
Specificity	
Immunogen Sequence	