

CCAR1 protein (STJ505500)

STJ505500

GENERAL INFORMATION

Product Type	Proteins
Short Description	CCAR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55749
Gene Symbol	CCAR1
Uniprot ID	CCAR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 645-695 on human Cell division cycle and apoptosis regulator protein 1.
Immunogen Region	645-695
Specificity	
Immunogen Sequence	