

DYRK1 protein (STJ505502)

STJ505502

GENERAL INFORMATION

Product Type	Proteins
Short Description	DYRK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1859
Gene Symbol	DYRK1A
Uniprot ID	DYR1A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Dual specificity tyrosine-phosphorylation-regulated kinase 1A.
Immunogen Region	300-350
Specificity	
Immunogen	
Sequence	