

DYRK2 protein (STJ505505)
STJ505505**GENERAL INFORMATION**

Product Type	Proteins
Short Description	DYRK2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8445
Gene Symbol	DYRK2
Uniprot ID	DYRK2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	