

5HT2A Receptor protein (STJ505508)

STJ505508

GENERAL INFORMATION

Product Type	Proteins
Short Description	5HT2A Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3356
Gene Symbol	HTR2A
Uniprot ID	5HT2A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human 5-hydroxytryptamine receptor 2A.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	