

5HT2C Receptor protein (STJ505510)

STJ505510

GENERAL INFORMATION

Product Type	Proteins
Short Description	5HT2C Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3358
Gene Symbol	HTR2C
Uniprot ID	5HT2C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-458 on human 5-hydroxytryptamine receptor 2C.
Immunogen Region	400-458
Specificity	
Immunogen Sequence	