

ARMET protein (STJ505512)

STJ505512

GENERAL INFORMATION

Product Type	Proteins
Short Description	ARMET-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
Uniprot ID	MANF_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 80-120 on human Mesencephalic astrocyte-derived neurotrophic factor.
Immunogen Region	80-120
Specificity	
Immunogen Sequence	