

IKAP protein (STJ505518)
STJ505518**GENERAL INFORMATION**

Product Type	Proteins
Short Description	IKAP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	230233
Gene Symbol	Elp1
Uniprot ID	ELP1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1280-1330 on mouse Elongator Complex Protein 1.
Immunogen Region	1280-1330
Specificity	
Immunogen Sequence	