

FOLR1 protein (STJ505521)
STJ505521**GENERAL INFORMATION**

Product Type	Proteins
Short Description	FOLR1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14275
Gene Symbol	Folr1
Uniprot ID	FOLR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	