

Apelin Receptor protein (STJ505526)

STJ505526

GENERAL INFORMATION

Product Type	Proteins
Short Description	Apelin Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	187
Gene Symbol	APLNR
Uniprot ID	APJ_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-377 rat Apelin receptor protein.
Immunogen Region	330-377
Specificity	
Immunogen Sequence	