

D4S234E protein (STJ505538)

STJ505538

GENERAL INFORMATION

Product Type	Proteins
Short Description	D4S234E-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25247
Gene Symbol	Nsg1
Uniprot ID	NSG1_RAT
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human Neuron-specific protein family member 1.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	