

Von Hippel Lindau protein (STJ505541)

STJ505541

GENERAL INFORMATION

Product Type	Proteins
Short Description	Von Hippel Lindau-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7428
Gene Symbol	VHL
Uniprot ID	VHL_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	
Immunogen Sequence	