

NSG2 protein (STJ505542)

STJ505542

GENERAL INFORMATION

Product Type	Proteins
Short Description	NSG2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51617
Gene Symbol	NSG2
Uniprot ID	NSG2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 130-171 on human Neuron-specific protein family member 2.
Immunogen Region	130-171
Specificity	
Immunogen Sequence	