

NSG2 protein (STJ505543)
STJ505543**GENERAL INFORMATION**

Product Type	Proteins
Short Description	NSG2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51617
Gene Symbol	NSG2
Uniprot ID	NSG2_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	
Immunogen Sequence	